

## Bariatric Battery-Powered Lift

HCPCS: E0635



- To charge, simply plug in unit Battery does not need to be removed to be charged
- Lift will not operate while plugged in Ideal for home-care use
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Base legs adjust easily and lock securely into open position with padded shifter handle
- Extra wide coated cradle with 360° rotation with 4-point hook up (Can fit under bed)
- Weight Capacity 600 lbs

PRODUCTS	Description
13244	Bariatric Lift,Elect,Homecare Style,1/cs

SPECIFICATIONS	Description
Base Clearance	2"
Base Height	4.5"
Base Width	45.75" (Open)   25" (Closed)
Boom Height	32" - 75.5"
Carton Shipping Weight	148 lbs.
Casters	3" Front   4" Rear
Motor	Linak Motor
Warranty	Lifetime on frame   1 year on components and electrical parts
Weight	132 lbs.
Weight Capacity	600 lbs

PARTS	Description
CPCFUSE	Fuse,10A,250V,CPC,13240,1/ea, EA
13244BMB	Base to Mast Bolt, 13244, 1/ea, EA
13240HC	Hand Control for 13240,13244,13246,1/ea, EA
13244F	Footplate for 13244,1/ea, EA
13240A	Actuator for 13240,1/ea, EA
13244MHB	Mast Hinge Bolt,13244/13245,1/ea, EA
13053PB	CradlePivotBolt13244/13245,1ea, EA
13244PH	Push Handle for 13244,1/ea, EA
13240BC	Battery Cell only for 13240,1/ea, EA
13240FP	Square Frame Plug,13240,1/ea, EA
13244SB	SpreaderBar for Base,13244,1/ea, EA
13240SBH	Handgrip,SpreaderBar,13240,1/ea, EA
13244BAR	Base Actuating Rod,13244,1/ea, EA

[Click here to view detailed product information and pricing.](#)

**Call 800-397-5899 to order today!**

**PARTS****Description**

13244CAL	Caster w/Lock,13244,1/ea, EA
13244ML	Lower Mast,13244,1/ea, EA
13244CPL	Caster Position Lock,13244,1/ea, EA
13240CA	Caster w/o Lock,13240,1/ea, EA
13244B	Base Only,13244 Lift,1/ea, EA
13240PC	Power Cord for 13240,13244,1/ea, EA
13244RPA	Receiver Plate Assm.13244,1/ea, EA