

## Sit To Stand Lift

HCPCS: E0635



- 450 lbs. Weight Capacity
- Dual controls - never have a down lift because of inoperable hand pendant
- Audible and Visual Alert if exceeds weight capacity of 450 lbs
- Object Obstruction Sensor - lift stops if detects patient or caregiver is in the way
- Adjustable knee pad
- Removable foot plate
- Lowest floor clearance in class
- Vertical and horizontal position ergonomic handles
- Battery can be charged on and off the lift
- Audible/Visual alert when battery is running low
- Only STSP450 power base includes
  - LCD screen for visual ready out
  - Maintenance monitoring
  - Tracking to support no lift policies
- BEST in Class Warranty: Frame Limited Lifetime,
- Electrical 5 Years, Battery 1 Year

### PRODUCTS

### Description

STSM450	Sit To Stand - Manual Base
STSP450	Sit To Stand - Power Base

### SPECIFICATIONS

### Description

Waranty	Limited Lifetime (Frame), 5 Years (Electrical), 1 Year (Battery)
Weight Capacity	450 lbs.

[Click here to view detailed product information and pricing.](#)

**Call 800-397-5899 to order today!**