

ROBOOTER

R E B O O T Y O U R L I F E



A SMART WHEELCHAIR WITH
THE ULTIMATE EXPERIENCE

www.roobooter.com



www.roobooter.com

ROBOOTER E40

ULTIMATE EXPERIENCE
NUMEROUS INGENUITY
ALL IN ROBOOTER E40



ENHANCED PERFORMANCE POWERFUL & AFFORDABLE

The cruising mileage 23km with 20ah battery, full driving force

APPEARANCE BEYOND DESCRIPTION ALL IN ROBOOTER E40



Ergonomic

Leather stitched armrest

High-end squeeze casting process, wear well

Fine matte spray treatment technology

COOPERATED SOFTWARE AND HARDWARE BEYOND YOUR IMAGE

Cruise mode, drive easily
Upgrade the system by app
Fault self-check and protection



CLASSIC WHITE



CLASSIC GREEN

PRODUCT PARAMETERS

Unfolded Size(Mm)	1000 (L)X 624 (W) X 930 (H) MM	Barrier Crossing Height	4CM
Folded Size(Mm)	700 (L)X 624 (W) X 450 (H) MM	Max Ditch Width	10CM
Seat Size(Mm)	700 (H)X 440 (W) X 440 (D) MM	Min Turning Radius	450MM
Cushion Size(Mm)	440 (L) X 440 (W) MM	Gradeability	9°
Folding Method	MANUALLY FOLDING	Front Wheel Material	HONEYCOMB PU TIRE
Max Weight Capacity	150KG	Rear Wheel Material	PNEUMATIC TIRES
Battery Capacity	20AH	Product Net Weight (Kg)	27KG
Endurance Mileage	23KM	Battery Weight	3.4KG
Max Speed	8KM/H	Color	WHITE、GREEN