

Innovation for people on the go

DreamStation Go provides the portability and packability travelers want and the effective therapy they need.



Charge

Onboard USB port allows for easy charging of mobile devices.



Power

Built-in power supply reduces the number of parts for you to manage for added convenience.



Safe

Complies with FAA requirements for in-flight use.¹



Practical

Bluetooth is standard with every device.



You're In control

Intuitive color touchscreen offers quick, simple setup and easy navigation.



Compatibility

Device is compatible with our DreamMapper monitoring self-management tool.



¹RTCA/DO-160G section 21, category M.



© 2017 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Caution: US federal law restricts these devices to sale by or on the order of a physician.

www.philips.com/respironics

RRDPGH MK 3/14/17 MCI 4107896 PN 1134215
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 · 724-387-4000

PHILIPS

DreamStation Go



Anywhere you go

DreamStation Go

Packs performance anywhere you go



DreamStation Go

Whether you're traveling for business or venturing out on the vacation of a lifetime, DreamStation Go provides an easy, reliable and portable PAP experience when you refuse to compromise.

Our sleek design has 1/3 fewer components than other travel PAPs,¹ and at half the size of our trusted System One and DreamStation platforms, DreamStation Go offers the same therapy and comfort features enjoyed by over 5 million PAP users.

¹When compared to against HDM Z1 and Somnetics Transcend devices (2016)

Convenience without compromise

Wherever your travels take you, sleep apnea follows. That's why DreamStation Go is an indispensable travel partner, packing powerful performance into a light, streamlined package. It even has a convenient USB port to charge your mobile device.

- Our smallest and lightest tubing¹
- Our smallest PAP machine
- Our lightest PAP at less than two pounds
- Complies with FAA requirements²

¹When compared to Philips 22mm performance tubing
²RTCA/DO-160G section 21, category M.



CPAP therapy even 'off the grid'

Overnight battery

When traveling to locations off the power grid, it can be difficult for those with sleep apnea to maintain their therapy. That's why we developed the DreamStation Go overnight battery.



DreamStation Go set at 10cm of pressure displayed an average run time of 13 hours¹ when running off the DreamStation Go overnight battery.

¹PAP run-time test conditions: CPAP mode, pressure 10 cm H2O, 12mm tubing, 37LPM of leak, 73.4 F room at an elevation of 1650 ft.



Small travel kit

Kit includes an easy-to-pack protective case for the DreamStation Go and a separate, washable bag for your mask and tube.



Medium travel kit

Kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes separate, washable bag for your mask and tube.

Learn more

Visit sleepapnea.com/dreamstationgo or call 1-844-240-1649.