

Aluminum Transport with 12" rear wheels

Features

- 12" (30.5 cm) wheels and handbrakes on top
- Powder coated aluminum frame is lightweight and durable
- Back folds down for easy storage and transport
- Seat belt for safety
- Comfortable nylon upholstery
- Full length permanent armrest
- Detachable footrests
- Larger rear wheels for better performance on uneven outdoor surfaces
- Weighs only 22 lbs (10 kg)!

Accessories

| Item Number | Description |
|--------------|-----------------------------|
| MDS85183FT* | I.V. Pole |
| MDS85181FT* | Oxygen Holder |
| MDS85190FT* | Oxygen / I.V. Combo |
| MDS85189FT | Anti-Tip Device |
| MDS85196 | Anti-Fold Anti-Theft Device |
| WCA806991 | Tinnerman Legrest Locks |
| MDSCHAIRCASE | Carrying Case with Strap |

*Recommended for use with Anti-Tip Device

300 lb | 136 kg
weight capacity >



MDS808210ARE

-  Red
MDS808210ARE
-  Blue
MDS808210ABE

| Item Number | Color | Seat W x D: Inch (cm) | HCPCS |
|--------------|-------|-----------------------|-------|
| MDS808210ARE | Red | 19 x 16 (48 x 41) | E1038 |
| MDS808210ABE | Blue | 19 x 16 (48 x 41) | E1038 |

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!