

MICROFLEX®

# BECAUSE SIZE MATTERS


The first ambidextrous  
exam glove available  
in ½ sizes.

Wearing exam gloves  
that are too loose or  
too tight may lead to  
hand fatigue and possibly  
career impacting injuries.

Your hands are special.  
Your gloves should be too.

**Microflex® Ultraform®** comes in a  
range of ½ sizes for a perfect fit,  
because no two hands are exactly  
the same.

If your glove isn't **Microflex Ultraform**,  
then it isn't half as good.



PLACE YOUR HAND HERE TO  
APPROXIMATE YOUR SIZE.

REQUEST A SAMPLE TO KNOW  
YOUR EXACT SIZE.

- Extra-Small 77 mm
- Extra-Small/Small 82 mm
- Small 88 mm
- Small/Medium 92 mm
- Medium 97 mm
- Medium/Large 102 mm
- Large 106 mm
- Extra-Large 115 mm

**NEW!**

**NEW!**

**NEW!**



Microflex®  
**Ultraform®**


**THE BEST FIT** JUST GOT BETTER  
Visit [Microflex.com](http://Microflex.com) for a sample

**Ansell**

## Microflex® Ultraform® Powder-Free Nitrile Examination Gloves

Microflex Ultraform with ERGOFORM technology, an innovation that reduces stress on joints, ligaments and tendons during repetitive tasks. Users experience more comfort and complete tasks with less fatigue minimizing the risk factors that lead to workplace injuries.

- Amazing fit and feel
- Now available in ½ sizes
- Ergonomic Certification
- Bulk Packed 300 gloves/box\*  
\*XL 250 gloves/box



## ERGOFORM™ Ergonomic Design Technology



PRODUCT INFORMATION		
Material	Powder-Free Nitrile	
Sterility	Non-Sterile	
Usage	Single Use Only	
Color	Cobalt	
Cuff Length	Standard 9.5"	
Finish	Textured Finger Tips	
ORDERING INFORMATION		
Size	Part #	Packaging
XS	UF-524-XS	300 gloves per box
XS/S	UF-524-XS/S <b>NEW!</b>	
S	UF-524-S	
S/M	UF-524-S/M <b>NEW!</b>	10 boxes per carton
M	UF-524-M	
M/L	UF-524-M/L <b>NEW!</b>	250 gloves per box
L	UF-524-L	
XL	UF-524-XL	

Visit [Microflex.com](http://Microflex.com) to learn more or request FREE samples